



12029-1

BARIATRIC
600 LB
WEIGHT CAPACITY

Deluxe Bariatric Bath Bench with Cross Frame Brace

Available
1 Pack

- 12029-1 With Tool Free Removable Back, Fully Assembled, 1/cs
- 12029-2 With Tool Free Removable Back, Fully Assembled, 2/cs
- 12028T Replacement Tips for 12028, 12029, 4/bx

FEATURES:

- 600 lb. Weight Capacity.
- Frame has width/depth cross braces attached with aircraft-type rivets.
- Large 16" seat depth.
- Durable, blow molded plastic bench.
- Tool free removable back.
- Aluminum frame is lightweight, durable and corrosion proof.
- Drainage holes reduce slipping.
- Limited Lifetime Warranty.



ASSEMBLED WITH TOOL FREE REMOVABLE BACK



EURO-STYLE CLIP FOR HEIGHT ADJUSTMENT!



12021-4

BARIATRIC
500 LB
WEIGHT CAPACITY

Deluxe Bariatric Bath Bench

Available
1 Pack

- 12021-4 With Back, Assembled, 4/cs
- 12021KD-1 With Back, Knock Down, 1/cs
- 12022-4 Without Back, Assembled, 4/cs
- 12022KD-1 Without Back, Knock Down, 1/cs
- 12021C Replacement Euro-Style Clip, 1/ea
- 12021T Replacement Tips for 12021, 4/bx

FEATURES:

- 500 lb. Weight Capacity.
- Frame has width/depth cross braces attached with aircraft-type rivets.
- Durable, blow-molded plastic bench.
- Aluminum frame is lightweight, durable and corrosion proof.
- Drainage holes reduce slipping.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12029-1/2	22"	17.25"	18.75"	16"	20"	16"	10 lbs.	11/21.5 lbs.	600 lbs.
12021-4	18"	20"	17.5"	12"	19.5"	14"	7.5 lbs.	32.5 lbs.	500 lbs.
12021KD-1	18"	20"	17.5"	12"	19.5"	14"	7.5 lbs.	8.5 lbs.	500 lbs.
12022-4	18"	20"	17.5"	12"	19.5"	14"	6.0 lbs.	27 lbs.	500 lbs.
12022KD-1	18"	20"	17.5"	12"	19.5"	14"	6.0 lbs.	8.0 lbs.	500 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.