



Ankle-Foot Orthosis 28U9:

- minimal weight
- excellent hygienic properties
- high Quality
- additional calf- and shin bone pads



Ankle-Foot Orthosis Comfort with a longer sole

This AFO from Otto Bock aids in swing phase toe clearance and reduces plantar flexion impact in early stance. Its thin-walled polypropylene construction has been optimized for increased resilience. The full sole plate provides a more precise guidance for the entire foot as well as enhanced pressure distribution. Thermoforming and machine grinding may be used to custom fit the orthosis to each individual. The 28U9 is notable for its minimal weight, high quality, and effective function.

Available for shoe sizes 35-37, 37-39, 39-41, and 41-44.



Ankle-Foot Orthosis 28U9

Indication

- Weakness of dorsiflexion
- peroneal paralysis



Fax order form

28U9 Ankle-Foot Orthosis -longsole-, Material: polypropylene, injection molded

Art.-No.	side	for shoesize	height	Quantity
28U9=L35-37	left	35-37	300 mm	____pc.
28U9=R35-37	right	35-37	300 mm	____pc.
28U9=L37-39	left	37-39	330 mm	____pc.
28U9=R37-39	right	37-39	330 mm	____pc.
28U9=L39-41	left	39-41	345 mm	____pc.
28U9=R39-41	right	39-41	345 mm	____pc.
28U9=L41-44	left	41-44	390 mm	____pc.
28U9=R41-44	right	41-44	390 mm	____pc.

Commission: _____

Customer-No.: _____

Sender/Receiver of goods

Otto Bock[®]

QUALITY FOR LIFE