



# Sensa Line

With SKINGUARD® Technology

- ⊕ Antibacterial
- ⊕ Neutral Odour
- ⊕ Skin-Friendly
- ⊕ Climate-regulating
- ⊕ Pleasantly soft
- ⊕ Latex-free

The Sensa Line active supports has now been equipped with SKINGUARD® technology. It contains a new kind of antibacterial fibre that prevents reproduction of bacteria and thus efficiently reduces odours. This results in significantly improved hygiene.

*Otto Bock*®

QUALITY FOR LIFE

# Size Selection


Sensa Line with SKINGUARD® Technology



## 50S5 Malleo Sensa

Side	Size	Circumference above the ankle (cm)
left/right	XS	17 - 19
left/right	S	19 - 21
left/right	M	21 - 23
left/right	L	23 - 25
left/right	XL	25 - 27
left/right	XXL	27 - 29

Colours: 0= beige; 1=grey



## 50K15 Genu Sensa

Size	D circumf. (cm)	F circumf. (cm)	F V-MAX circumf. (cm)
XXS	25 - 28	35 - 38	-
XS	28 - 31	38 - 41	-
S	31 - 34	41 - 44	-
M	34 - 37	44 - 47	47 - 50
L	37 - 40	47 - 50	50 - 53
XL	40 - 43	50 - 53	53 - 56
XXL	43 - 46	53 - 56	56 - 59


Colours: 0= beige; 1=grey  
Can be worn either on the left or on the right knee



## 50R10 Lumbo Sensa 50R11 Lumbo Sensa, women

Size	Body circumference in cm (approx. 2 cm above the iliac crest)
XS	70 - 80
S	80 - 90
M	90 - 100
L	100 - 110
XL	110 - 120
XXL	120 - 130

Colours: 0=beige; 1=grey  
50R10 Lumbo Sensa (straight style)  
50R11 Lumbo Sensa, women (waisted style)



## 50A7 Epi Sensa

Size	Forearm circumference (cm)
XXS	17 - 19
XS	19 - 21
S	21 - 23
M	23 - 25
L	25 - 27
XL	27 - 29
XXL	29 - 31

Colours: 0=beige; 1=grey  
Can be worn either on the left or on the right elbow



## 50P13 Manu Sensa

Side	Size	Wrist circumference (cm)
left/right	XS	14 - 15
left/right	S	15 - 16
left/right	M	16 - 17
left/right	L	17 - 18
left/right	XL	18 - 19
left/right	XXL	19 - 20

Colours: 0=beige; 1=grey

## Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt/Germany · Phone +49 5527 848-1424  
Fax +49 5527 72330 · export@ottobock.de · www.ottobock.com