



Otto Bock®

QUALITY FOR LIFE



## 50A10 Omo Immobil

Shoulder Abduction & External Rotation Control in One

The Omo Immobil gives you a single solution for patients who need precise shoulder positioning due to multiple surgical and non-surgical indications. The Omo Immobil makes it simple to provide abduction and external rotation so you have the options you need in a single Orthosis. In addition to its wide array of applications, the Omo Immobil's comfortable design helps increase patient compliance.

USA 800.328.4058 • Fax 800.962.2549 • [www.ottobockus.com](http://www.ottobockus.com)



# 50A10 Omo Immobil

## Shoulder Abduction & External Rotation Control in One

©2008 Otto Bock HealthCare LP • 08091538.1 • 11/08

The 50A10 Omo Immobil shoulder orthosis features numerous options for post-operative and post-traumatic shoulder indications and includes multi-directional shoulder positioning.

### Features and Benefits

- Easily adapts for left or right shoulder, which helps you reduce stock volume
- Five abduction angles options: 0°, 15°, 30°, 60° or 90°
- Three external rotation options: 10°, 20° or 30°
- Simplified patient fittings because there are no tools required
- Low profile, slim design
- Breathable, washable, removable cover for increased patient compliance
- Waist belt reduces load of the shoulder strap on the neck for greater patient compliance
- External rotation can be used in any combination with all abduction angles

### Patient Indications – Omo Immobil 50A10

- Rotator cuff surgery
- Non-operative treatment of anterior shoulder dislocation
- Subacromial decompression
- Total shoulder arthroplasty
- Bankhart repairs

### Patient Indications – Omo Immobil External Rotation Set 29A222

All indications requiring immobilization of the shoulder at 10°, 20° or 30° of external rotation

### Sizing Chart

The size of the orthosis is selected according to body height.

Material No.	Size	Body Height
50A10=	S	150 – 165 cm / 58.8 – 64.5 in
50A10=	M	165 – 175 cm / 64.5 – 68.5 in
50A10=	L	175 cm / 68.5 in
29A222	S, M, L	External Rotation Set



Omo Immobil



Omo Immobil Rotation

### Ordering Information

For more information about the Omo Immobil or Omo Immobil Rotation, talk with your local Otto Bock Sales Representative at 800.328.4058 (U.S.) or 800.665.3327 (Canada). You can also find more information online at [www.ottobockus.com](http://www.ottobockus.com) or [www.ottobock.ca](http://www.ottobock.ca).

*Otto Bock*<sup>®</sup>

QUALITY FOR LIFE