

Elite Cushion

Future Mobility Healthcare manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Renown for quality and service we bring innovation to all of our products for the home healthcare market


The Prism Elite Cushion was designed for users that require moderate postural support and offers gentle contours helping to promote correct seating posture.

A top layer of soft, comfortable poly foam, along with the visco and high density polyurethane foam offers pressure distribution and relief. We truly have designed an Elite cushion.



MADE IN CANADA

Elite Cushion



CUSHION THICKNESS

- A—Height at adductor = 4"
- B—Height at abductor = 4 1/2"
- C—Height at rear IT = 3"
- D—Leg trough = 3"

SPOTLIGHT ON

- Waterfall front edge for foot propelling
- Lightweight cushion
- Non-slip bottom on washable incontinent cover
- Core layer of visco-elastic foam to provide comfort and pressure reduction

WARRANTY

Two (2) years limited on Cushion and 90 days on cover.

PRODUCT WEIGHT

2.5 lbs

USER WEIGHT

250 lbs

<i>Stock sizes—same day shipments</i>			
ADP CODE: SESFM0060			
16" x 16"	EC1616	16" x 18"	EC1618
18" x 16"	EC1816	18" x 18"	EC1818
20" x 16"	EC2016	20" x 18"	EC2018
<i>Non-stock sizes—shipments within 3 business days</i>			
ADP CODE: SESFM0060			
16" x 20"	EC1620	17" x 16"	EC1716
17" x 18"	EC1718	17" x 20"	EC1720
18" x 20"	EC1820	20" x 20"	EC2020
<i>Custom sizes—shipments within 7 business days</i>			
ADP CODE: SESFM0065			
ECCUS			



3223 Orlando Drive, Mississauga, ON, L4V 1C5, 905-671-1661, 888-737-4011, fax 866-670-3377 www.future-mobility.com
 Future Mobility Healthcare Inc. Orion and Prism are trademarks of Future Mobility Healthcare Inc.