

## Some DREAM Series Options

Colors in this document are approximations; we do not guarantee color accuracy from printed or electronic images.



Raised panel doors to reduce sensory input

Natural finish on oak

Interior padding  
Window in Head




Rosewood stain on maple



Drawers  
Access port

Michael's cherry stain on maple

Lathe turned finial posts



Accent Carvings and Engraving (below)



Fruitwood stain on maple

## Safety Mesh Canopy



When high sides aren't enough consider adding a safety top. This canopy is made from non stretchable mesh with marine grade coated nylon trim for strength. This durable fabric is stretched over an incredibly strong aluminum frame which is bolted to our corner posts.

When installed the canopy adds approximately 15" of additional height to the bed. The doors on the bed continue to operate freely even when the canopy is in place. The canopy itself does not open.

Choose from over a dozen trim colors.

The BBG canopy can be attached to any series or model of BBG bed, but needs to be ordered with the bed and currently cannot be retrofit.



**BBG Beds by George**®  
Safety Beds for the Medically Fragile



**BBG Beds by George**®



25416 CR 6 Ste 106, Elkhart, IN 46514  
Phone 877-286-1994  
Email: info@bedsbygeorge.com  
www.bedsbygeorge.com

Proudly made in the USA

Everyone sleeps better with  
a BBG bed in the house.

www.bedsbygeorge.com

Pediatric  
 Teen  
 Adult  
 Geriatric

# More Choices with 2 Series of Beds!

## All Beds Include

- Handcrafted Bed Frame
- Impact resistant door panels
- Mattress w/ water resistant, anti bacterial, low shear cover
- Bunkie board rated at 300 pounds
- Solid wooden construction, no plastic edging or trim.
- Massive solid wood corner posts for stability and strength
- Doors fully open along entire side and store out of the way
- Locking Doors
- Meets or exceeds FDA standards for entrapment and entanglement
- UL and CSA approved
- 5 year limited warranty



## Standard Door System

## Hi Side Door System

**Slumber Series**  
 Basic Safety



### Available with:


1. Non adjustable mattress, Fixed height.
2. Adjustable mattress (crank) , Fixed height.
3. High Low, adjustable mattress (electric)




**Dream Series**  
 Personalized Safety



**Reduced gap to prevent burrowing, entrapment or tunneling under mattress**



- Doors 22" above mattress
- Transfer height 25"
- Doors on front and back standard



- Doors 32" above mattress
- Transfer height 25"
- Dual front doors open together
- Back Door optional

Patent Pending

Add options such as canopy, interior padding, iv poles, wheels, access ports. (even more options in Dream Series!)