

Marquis 888^{SLN}

SHOPPRIDER

Mid-size scooter

Shoprider's Marquis 888^{SLN} is part of our classic line of scooters with the perfect blend of style and four-wheel stability performance.

- > Shoprider 4-wheel scooter
- > Full light package
- > Full suspension, anchor hooks
- > 50AH batteries with 4A charger
- > Convenient battery access panel
- > Delta Tiller
- > Available colours: red, black or pacific blue



888^{SLN}

Specs	Weight capacity	Dimensions LxWxH	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v chrgr	Driving dist.	Turning radius	Safe climb ^o
888 ^{SLN}	136KG/300LBS	47.5" x 22.5" x 42"	99KG/218LBS	0.87HP	8KPH/5MPH	2 x 50AH	F10" R10"	4AMP/HR	40KM/25M	1460MM/58"	8°/14%

All specifications subject to change without notice.

Dealer stamp



eclipse[™]
medical inc.

Discovering your world again—in style

www.eclipsemedical.com