

permobil



AN ORGANIC EXPERIENCE.

The C400 Corpus® combines innovative design, styling and advanced engineering concepts. The C400 maintains a small profile but combines better indoor maneuverability with superior outdoor performance. With a rugged, advanced suspension system, the C400 has advanced obstacle climbing capabilities. The Corpus seating system of the C400 is modular in design, allowing for growth and added seat functions. Ergonomically pleasing with power tilt and recline, customization of the Corpus is almost limitless...an important feature when you are becoming one with nature.

C400
CORPUS



POWER FUNCTIONS



Power adjustable tilt



Crystal Blue



Vintage Wine



Metallic Green



Xtreme Yellow



Black Onyx

SHROUD COLORS



R-Net electronics with ICS



Power adjustable recline



Power seat height adjustment



Power adjustable legrests

SPECIFICATIONS

Maximum Speed	4.5 MPH
Driving Range	19-25 Miles
Max User Weight	250 lbs
Ground Clearance	2.8"
Battery Capacity	Group 24
Width	24" - 28.5"
Length	34.5"
Recommended Battery	Group 24
Gradient Capacity	10°
Power Seat Elevator	8"

Seat Widths	17" - 23" (Adjustable)
Seat Depths	14" - 22" (Adjustable)
Backrest Heights	16" - 22"
Seat To Floor Height	19.25" - 23.5"
Electronics (standard)	R-Net
Alternative Drive Controls	Yes w/R-Net
Color Options	Black Onyx Metallic Green Xtreme Yellow



THE POWER OF MOBILITY

Permobil Inc.
 6961 Eastgate Boulevard
 Lebanon, Tennessee 37090
 800.736.0925
www.permobil.com

C400042008